

| Technical Guidance (PaX-i3D Guide for Defect Prevention of AnyGrabber Card 5) | | | | Creator | Confirmer | Approver |
|--|--|----------|------------------------|---------|-----------|----------|
| | | | | | | |
| Creator | Eunice Lyou | Doc. No. | GCS-TECH-20140114B-ENG | | | |
| To | | Cc | | | | |
| Object | PaX-i3D | Category | AnyGrabber Card | | | |
| Subject | PaX-i3D Guide for Defect Prevention of AnyGrabber Card 5 | | | | | |

| Revision History | | | |
|------------------|----------|----------------|--------------------|
| Date | Revision | Editor | Comments |
| 2014.01.14 | 0.0 | Augustine Jeon | Created |
| 2014.01.14 | 0.0 | Augustine Jeon | Conversion |
| 2014.01.22 | 0.0 | Eunice Lyou | Conversion Revised |

Contents

| | | |
|----------|--|----------|
| 1 | General Provision..... | 2 |
| 1.1 | Objective | 2 |
| 1.2 | Application Scope..... | 2 |
| 2 | Description | 2 |
| 2.1 | Symptom | 2 |
| 2.2 | Version Information | 2 |
| 3 | Work Process | 2 |
| 3.1 | Update of the Driver of AnyGrabber Card 5..... | 2 |
| 4 | Attachments..... | 5 |

1 General Provision

1.1 Objective

- This is to provide the guide for defect prevention of AnyGrabber Card 5 of PaX-i3D.

1.2 Application Scope

- Model : PaX-i3D
- S/N: Applicable for the product released before S/N 052-0875

2 Description

2.1 Symptom

- 1) PC reboots after the occurrence of the blue screen in PC screen when selecting the ready button in capture software.

2.2 Version Information


- 1) Version Information of Driver

| Classification | Before Modification | After Modification |
|----------------|---------------------|--------------------|
| AG50 Driver | V 1.0.3.4 | V 1.0.3.5 |

3 Work Process

3.1 Update of the Driver of AnyGrabber Card 5


- 1) Uninstall the AnyGrabber3 N1000(x64) at the Device Manger → Imaging devices tab.




- 2) Check on the check box of “Delete the driver software for this device” and press OK.




- 3) Multimedia Controller” is searched automatically at the ‘Other devices tab’ of Device Manager.




- 4) Move to Control Panel → Program → Programs and Features and uninstall the “AnyGrabber3(5)(x64)”.




- 5) After unzipping the "Frame_Grabber_AG5.zip" file, move to the Driver tab and select the system type according to the customer's system.




- ① In case of system with 32bits




- ② In case of system with 64bits





- 6) Proceed the installation followed by Install Shield Wizard.





System must be rebooted after installing the Driver.

7) Confirm the Image device installation and check the version of AnyGrabber Card 5.



4 Attachments

- 1) Frame_Grabber_AG5.zip
- 2) PaX-i3D AG50 Driver Issue List.xlsx