

NEW

Endopositioner

Better vision in the root canal



- Avoiding shadows in the root canal, because the light beam and the visual axis do not cross - for optimum vision
- Usable with opt-on loupe
- Simple and fast attachment due to magnetflex-technology
- Combinable with the wireless light system spot-on^{®cordfree} and the spot-on^{®nxt}



Contact:

info@orangedental.de
tel: +49 7351 474 990
www.orangedental.de

orangedental
premium innovations




NEW


Endopositioner

Better vision in the root canal

Due to the angle α of the light beam a shadow („blind spot“) is created in the root canal, which impairs the view.



The endopositioner avoids shadows in the root canal, because light beam and visual axis do not cross - for optimum vision.



Contact:

info@orangedental.de
tel: +49 7351 474 990
www.orangedental.de

orangedental
premium innovations

