


FUSSEN by orangedental

Intraoral scanner model S6500

Model S6500



- Large field of view 16 x 13 mm
- Complete jaw scan in less than 2 minutes
- High accuracy only 10 µm tolerance
- Super light and compact only 220 g
- 3 standard and 1 mini scan tip
- New ScanFlex software with many AI functions



- New, modern and stable cart with mini PC, touch screen and emergency power supply

FUSSEN by orangedental

Intraoral scanner

Model S6500

	S6500
Type	Cable
Scanprinciple	Triangulation
Light source	Structured light
Complete jaw scan (min)	< 2
Scanning speed (FPS)	20-25
Depth of field (mm)	0-15
Field of view (FOV)	16x13/12x9
Anti-fogging (min)	< 1
Weight (g)	220
Size (mm)	240x45x36
G-Motion	Yes
Precision single crown	10µm (±1,1µm)
Size of the scanning tips (mm)	Standard: 95,3 x22,5x16,6 Mini: 94x 17,9x 13,2
Included scan tips	3 standard scanning tips 1 Mini scanning tip
RRP plus VAT.	12.900 €

